

Download/installation instructions:



1) To download the software yEd (freeware), follow the link:


<https://www.yworks.com/products/yed/download>

2) To download the decision tree file for the yEd software, open


https://azcollaboration.sharepoint.com/sites/AB866/_layouts/15/start.aspx#/

yEd/decision tree user instructions

| | |
|---|---|
|  <pre>graph TD; A{Formulation type} --> B[ER]; A --> C[IR]; A --> D[DR];</pre> | <ul style="list-style-type: none">• Overview of the decision tree• The color code indicates the three main branches:<ul style="list-style-type: none">○ ER (yellow)○ IR (blue)○ DR (green) |
|  | <ul style="list-style-type: none">• Clicking on the small "+" on the upper left-hand side reveals the next level of the decision tree• TO return to a higher level of the decision tree, click on the "-" |




- You can change the selected window frame by using the classical scrolling bars or by holding down the right mouse button and moving the mouse




- You can zoom in or out by either pressing on the magnifying glass with the "+" or "-", or by using the mouse wheel

- Sometimes a part of the decision tree which is not part of the opened sub-level overlays some parts of the decision tree (in this case the IR box). To only show the decision tree part of interest, single left click somewhere on the grey area of the more detailed decision tree part.



- The deepest level of the decision tree are the overview pages with information on the selected apparatus/methods. Currently, these pages are linked in the decision tree but are located elsewhere
- All boxes which direct the user via a hyperlink to another document are colored in blue
- To open an overview page, right mouse click on a box which says "...click here for overview", then select "go to URL" (in the example it is written in German), or mark the box and press "F8" on the keyboard



- If you want to quit the software, even if you haven't changed anything, the software will ask if you want to save the changes. Please select "No" (also in German on the screen shot)

Glossary:

ER = extended release formulation
 IR = immediate release formulation
 DR = delayed release formulation
 ASD = artificial stomach-duodenum model
 D/S = dose solubility ratio
 FaSSIF = fasted state simulated intestinal fluid
 FaSSGF = fasted state simulated gastric fluid
 TM = transfer model